

[Home](#)

[News](#)

[Healthcare Zones](#)

[Latest issues](#)

[Links](#)

[Email bulletin](#)

[Advertise](#)

Bookmark this site

[Click here to bookmark hospitalhealthcare.com](#)

Related sites:

[Hospital IT Europe](#)
[Hospital Imaging and Radiology Europe](#)

Latest News

[Tweet](#)

[Like](#)

2

Sample Integrity with the NEO Plate 96 Friday 31st August 2012


LATEST NEWS

The monoblock detachable NEO Plate 96

As NEO plate is monoblock, the tubes cannot be moved. That means you can still have full positive identification without the need for a scanner however with the amazing "Craker Pen" tool a single microtube can be extracted from the monoblock NEO Plate 96.

Neo Plate 96, with its 96 detachable microtubes and a SBS footprint, is resistant to most organic solvents and provides sample storage in -135°C mechanical freezers. This plate is available in different volumes, including 70µL, 500µL and 1mL. Select the optimum size plate to make the best possible use of your available freezer space.

For further information on NEO Plate 96, please call at +33 (0) 1 34 30 42 08 or email at marketing@biosamplingsystems.com



RELATED ARTICLES

[IMC's drug-monitoring systems NEO Sample Management BioSampling Systems and Ramcon A/S Announce Distribution Partnership](#)
[New inline pump reduces downtime and improves accuracy](#)

[Why do we need accreditation?](#)

USER TOOLS

[Have your say - comment on this item](#)

[Email to friend](#)

[Print-friendly version](#)